




Collapsible MDI Holding Chamber:

An Innovative Alternative

LiteAire®'s unique dual-valved MDI holding chamber design delivers pop-up convenience and effective drug output at a fraction of the cost.

In most clinical settings, the LiteAire MDI holding chamber can reduce costs by replacing existing rigid plastic holding chambers or inefficient spacers with a cardboard alternative. This unique design allows the LiteAire to be reused by a patient over multiple doses and meets and often exceeds the performance of plastic holding chambers.

LiteAire cross section with the side panel removed.



Convenient, individually packaged, 25-count dispenser box

- Dual Valved
- Pop-Up to Use
- Portable/Stores Flat
- Labeled for One Weeks Use
- No Natural Latex
- Printed with Soy Ink

The LiteAire delivers true value when used in:

Emergency Department/Outpatient Clinic: For patients requiring the use of a holding chamber for immediate treatment.

Pulmonary Function Testing: Clean, disposable, ready to use in a truly single patient use environment; reduces the risks of patient cross-contamination; designed to accept all MDI's.

Short-Term Respiratory Infections: A cost effective treatment method that meets the short-term needs of the acute respiratory patient.

Patient Compliance: Complies with single patient use labeling and regulations, enhanced infection control and greater patient compliance with added convenience and portability.

Ordering Information

Single 25-Count Dispenser Box (Part No TH1304)

Case of Four 25-Count Dispenser Boxes (Part No TH1306)

For more information, telephone Niche Medical on 1300 136 855 or email us at info@nichemedical.com.au

NicheMedical

Telephone: 1300 136 855 Facsimile: 1300 664 204
Email: info@nichemedical.com.au
www.nichemedical.com.au

